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Operating instructions for Orgone Zapper Standard, Orgonezapper El Silverado and Clarkzapper

Application:

The zapper is designed as a wellness device to relax muscle tension. It has not been tested or approved for healing illnesses. With the help of the elastic band, the Zapper can be worn anywhere on the arms or legs. It can also be applied to any other part of the body without the band. The electrodes should be in light contact with the skin.

How it works:

The zapper produces a weak "tingling current" of 9V pulsed direct current. This is found by many users to be relaxing and invigorating.

Safety:

The zapper works with very low currents of 20mA maximum, which excludes any danger. However, it should not be used by people with pacemakers or during pregnancy.

Battery replacement and change:

The Zapper is supplied without a battery for flight safety reasons. A 9V battery of type PP3 is required. To insert the battery, loosen the cover with your fingernail at the corresponding slots and lift it off. The battery must lie flat on the bottom of the unit without a cable underneath, because then the lid would not close. Likewise, the push button must be positioned on the side facing away from the battery. Make sure that the polarity of the battery is correct.

Switching on and off:

The red LED light flickers when the zapper is switched on and the battery is correctly installed.

Care and maintenance:

Discolouration of the skin electrodes is normal. Especially with the silver electrodes, oxidation takes place, which is known to you from silver cutlery as "tarnishing". Since it is pulsed direct current, this will always be more pronounced on one of the two electrodes. There is no need to polish the electrodes bright all the time. Preferably, you should only wipe the electrodes with diluted lemon juice with a damp cloth. Metal polishes are not necessarily skin-friendly and should only be used in extreme cases. Then remove all residues with a damp cloth. Please be gentle with the small slide switch and the zapper will serve you faithfully for many years.

Warranty:

If used properly, the usual warranty period according to the German Civil Code (BGB) applies. The manufacturer is not liable for improper use other than described, negligent or intentional misuse or modification of the device.

Disposal:

The Zapper is considered electronic waste and must not be disposed of with household waste. It must be taken to a suitable waste disposal point. Batteries must be taken to battery collection points. The paper portion of the packaging should be sent to waste paper recycling.

Manufacturer:

Orgone Power (Pty) Ltd
166 Frances Street
Observatory 2198
ZA Gauteng
Republic of South Africa
orgonisehelpdesk@gmail.com
Tel +27 102 370180

Importer and agent for the European Union:

Mario Bürge
Neufahrwasserweg 2
DE 12685 Berlin
Germany
buergemario@gmail.com
Tel +49 30 92215822